



NEW CLOE is a recliner with a refined and elegant design, ideal for intensive use in lounge and relaxation areas. It features a reinforced solid wood internal frame upholstered with high-density polyurethane foam.

It is the perfect addition to any lounge or lobby setting. The seat and footrest feature continuous upholstery, enhancing safety during use.

Ideal for creating warm, homelike environments.

TECHNICAL DATA SHEET

NEW CLOE

Design: **Seniorcare Estudio**



TECHNICAL DATA SHEET

NEW CLOE

Design: Seniorcare Estudio

RECOMMENDED USE

Suitable for lounge areas and bedrooms.

SEATS AND BACKRESTS

The internal structure of the armchair is made of solid wood, with MDF side panels. The seat and backrest suspension system is based on high-durability NEA elastic webbing. Its high rubber content provides great resistance and excellent performance.

The entire unit is upholstered with high-density flexible polyurethane foam: 30 kg/m³ in the seats and 28 kg/m³ in the backrests, with HR resilience. The seat and footrest feature continuous upholstery to prevent gaps that could cause injuries.

ARMRESTS

The armrests feature an internal solid wood structure with particleboard side panels. They are upholstered with high-density flexible polyurethane foam (28 kg/m³) with HR resilience. In addition, removable upholstered armrests can be incorporated to facilitate the user's lateral transfer.

LEGS

Made of solid wood, available in natural, walnut or oak finishes.

MECHANISMS

Manual recline mechanism allows the user to manually adjust the coordinated inclination of the backrest and footrest.

Manual recline mechanism - GRAVITY ZERO position: allows the user to manually adjust the coordinated inclination of the backrest and footrest.

It includes the option to adopt the GRAVITY ZERO position, in which the legs are at the same level as or higher than the heart to promote blood circulation. This relaxation position helps relieve pressure on the lumbar and cervical muscles.

Motorized recline mechanism allows the user to adjust the desired inclination using a remote control. Powered by an electrical supply. By holding the button, the user can complete the full motion to the final reclined position or stop at any intermediate point. Once closed, the chair returns to the sitting position.

Motorized recline + lift mechanism allows the user to adjust the desired inclination using a remote control. Powered by an electrical supply.

By holding the button, the user can complete the full motion to the final reclined position or stop at any intermediate point. Once closed, the user can continue pressing the same button to complete the movement to the standing position.

*In the motorized versions, the optional GRAVITY ZERO position can also be included. In this position, the legs are at the same level or higher as the heart to promote blood circulation. This is a relaxation position that helps relieve pressure on the lumbar and cervical muscles.

TRANSPORT

The TRANSPORT version includes a transport handle and 75 mm diameter Agila wheels from Tente, of the highest quality, two of which are equipped with brakes, allowing the chair to be locked when not in motion.

The Nepal castors comply with the **UNE-EN 12530:1999 standard** –

Castors and wheels. Castors and wheels for manually propelled equipment used in public institutions.

TESTS PERFORMED. COMPLIANCE WITH THE STANDARD

Tests Performed at: **AIDIMME - Technological Institute of Metalworking, Furniture, Wood, Packaging and Related Industries**

The tests carried out on the model are those indicated and applicable according to the following standards:

- **UNE-EN 16139:2013 vc2015, Furniture - Strength, durability and safety - Requirements for non-domestic seating for general use under Level 1 public use conditions.**

The test procedures are described in the following standards:

- **UNE-EN 1728:2013, Furniture - Seating - Test methods for the determination of strength and durability,**
- **UNE-EN 1022:2019, Furniture - Seating - Determination of stability.**

The tests were carried out following the procedures described in the above standards, and the result was **COMPLIANT**.

Additionally, a durability test of the recline mechanism was performed according to **UNE-EN 13759:2012-1, Furniture - Mechanisms for the operation of seating and sofa beds**. After applying the prescribed number of cycles, no loss of functionality was observed.

ECOLOGY

ISO 14001 Certification

The production processes of SENIORCARE INTERIORS are certified in accordance with the ISO 14001 environmental standard, with the aim of minimizing any impact on the environment.

Seniorcare Interiors holds the Certificate of Registration in the Carbon Footprint, Offset and CO₂ Absorption Projects Registry of the Spanish Ministry for the Ecological Transition and the Demographic Challenge.

CERTIFICATES AND STANDARDS

Certificado ISO 9001

SENIORCARE INTERIORS tiene certificada su gestión de calidad conforme a la normativa ISO 9001 en todos los procesos de diseño, fabricación y comercialización de sus productos.

WARRANTY

5 years.

FINISHES

Wood: legs available in natural, walnut, or oak finishes.

Upholstery: all fabrics listed in our finishes book.

All our fabrics—except Apollo—are vinyl and feature antibacterial, antimicrobial, and antifungal properties.

They offer high resistance to abrasion, reaching up to 300,000 Martindale cycles in some cases.

WEIGHT CAPACITY

140 - 150 kilos

Información técnica / Technical information / Information technique



**Sillón relax CLOE asiento
continuo**

CL00011

H 108 x W 74 x D 95/170 cm
SH 46 cm | AH 62 cm



**Sillón relax CLOE
posición ZERO**

CL00020

H 108 x W 74 x D 95/170 cm
SH 46 cm | AH 62 cm



**Sillón relax CLOE motor +
reincorp. + posición ZERO**

CL00030

H 108 x W 74 x D 95/170 cm
SH 46 cm | AH 62 cm



**Aumento opción sill
traslado**

CLO00110



senior**care**

Camino de la Estación s/n
Polig. i-4, Parcela 5
03330 CREVILLENTE (Alicante)
Tels.: +34 96 540 53 30 / 96 540 35 12
Fax: +34 96 668 05 62
e-mail: info@seniorcare.es

www.seniorcare.es